



# Daniels

## Containment Systems

With clinically engineered containment systems for all healthcare wastestreams that ascribe to the highest levels of safety, security and environmental responsibility, Daniels Health is the industry leader in product & service that protects healthcare workers in the handling of high-risk waste.



### RMW SOFT & SHARPS WASTE

Clinically designed reusable containers for the collection of RMW and sharps waste



### CHEMOTHERAPY TRACE & BULK

Clinically designed reusable containers for the collection of trace chemotherapy waste Bulk chemotherapy container options available.



### HAZARDOUS WASTE

Hinged-lid containers that are leak proof and tamper proof.



### PHARMACEUTICAL WASTE & RCRA

Clinically designed reusable containers for disposal of drugs and Expired pharmaceuticals. Permanent security locks to prevent tampering and unauthorized access.



### SURGICAL DEVICES & INSTRUMENTS

We provide quarterly reports of percentage of waste diverted from landfill and the environmental impact.

# Sharpsmart Sharps Containers

## REGULAR



### 14-SERIES CONTAINER

**Shown:** Sharpsmart

2 Gallon Access Plus  
reusable Daniels collector

**Fill Capacity:** 8 Quarts

H 11inch x W 15.5 x D 8.3



### 22-SERIES CONTAINER

**Shown:** Chemosmart

4 Gallon Access Plus  
reusable Daniels collector

**Fill Capacity:** 18 Quarts

H 17.3 inch x W 15.5 x D 8.3



### 32-SERIES CONTAINER

**Shown:** Pharmasmart

8 Gallon Access Plus  
reusable Daniels collector

**Fill Capacity:** 28 Quarts

H 24 inch x W 15.5 x D 8.3



### 64-SERIES CONTAINER

**Shown:** Surgismart

16 Gallon Access Plus  
reusable Daniels collector

**Fill Capacity:** 56 Quarts

H 28.5inch x W 15.5 x D 8.3



### 28 GALLON TUB

**Shown:** RMW Tub

28 Gallon Daniels reusable tub

**Dimensions:**

H 22.9 inch x W 18.9 x D 18.9



### RCRA CONTAINER

**Shown:** Haz Waste Container

55lb Capacity Duo-port disposable  
container with hinged lid

**Dimensions:**

13 L x 17.3 W x 10.8 H