

The Daniels Utility Suite

100% CUSTOMIZABLE
100% POINT-OF-USE EFFICIENCY



Daniels' Utility Suite brings ultimate mobile efficiency to waste management and patient care. This clever interchangeable system of accessories enables multi-wastestream point-of-use disposal when wheeled into clinical environments or patient wards, as well as facilitating within-reach access to hand sanitizer, clean gloves, swabs, medication and small procedural items.

All Daniels utility accessories are compatible with the Daniels Accessmart trolley. Easy to assemble parts are supplied with screws that insert directly into the Daniels Accessmart's pre-engineered recesses, and dual side and rear mounting options allow optimal versatility in the customization of each unit.



Daniels Utility Accessories

EACH UNIT SOLD SEPARATELY



Accessmart Side Bracket | Enables dual mounting of an S14 or S22 container to the side of an Accesssmart trolley.

Product Code : 4001916000 | Dimensions : 4.6" (H) x 10.2" (W) x 1.2" (D)

Compatibility : S14, S22 Regular and Access Plus containers



Bag Grabber | Designed to hang a plastic bag to the side or back of an Accesssmart trolley, the Daniels bag grabber enables segregated waste disposal at the point of generation.

Product Code : 4006502000 | Dimensions : 3.0" (H) x 2.2" (W) x 1.1" (D)

Features : The bag grabber can be mounted on either the back, or the side of an Accesssmart trolley.



Small Glove Box Holder | Designed to mount glove boxes to the side of an Accesssmart trolley.

Product Code : 4006716000 | Dimensions : 14.0" (H) x 2.8" (W) x 3.1" (D)

Features : This product can be mounted to other surfaces as required. Holds up to 3 glove boxes, maximum dimensions are 2.6" (H) x 9.3" (W) x 4.8" (D).



Large Glove Box Holder | Designed to mount glove boxes to the side of an Accesssmart trolley.

Product Code : 4006616000 | Dimensions : 14.0" (H) x 2.8" (W) x 4.0" (D)

Features : This product can be mounted to other surfaces as required. Holds up to 3 glove boxes, maximum dimensions are 3" (H) x 9.3" (W) x 4.8" (D).



Hand Sanitizer Ring | Designed to firmly secure a bottle of hand sanitizer to the back of an Accesssmart trolley.

Product Code : 4006116000 | Dimensions : 2.5" (H) x 2.0" (W) x 2.0" (D)

Features : This product can be mounted to other surfaces as required. Compatible with most hand sanitizer brands. Maximum neck diameter of up to 1 inch.



Side Tray Kit | Designed to hold small procedural items for use in minor bedside treatments (i.e. kidney dish, clean syringe packaging and medication, gauze, bandages).

Product Code : 4004302000

Dimensions : 21.3" (H) x 12.5" (W) x 10.7" (D)

Features : Mounts to the side mounting point of an Accesssmart trolley.



Top Tray Kit | Designed to hold small procedural items for use in minor bedside treatments (i.e. kidney dish, clean syringe packaging and medication, gauze, bandages).

Product Code : 4004902000 | Dimensions : 23.2" (H) x 17.7" (W) x 14.0" (D)

Features : Mounts to both side mounting points of an Accesssmart trolley.



PICTURED : ACCESSMART TROLLEY WITH MEDISMART CONTAINER FRONT MOUNTED AND S14 SHARPSMART SIDE MOUNTED