

# Meet SECURE a SHARP

The solution for **secure, portable and convenient** sharps disposal. Designed for safety in healthcare settings.



## When to Use Secure a Sharp?

Secure a Sharp containers are designed for safe sharps disposal in non-acute healthcare settings. With a tethered cap, puncture-proof design, and secure locking mechanism, our containers are designed to protect healthcare workers and patients in the securing of used sharps.



### Step 1:

Dispose of used sharps into your Secure a Sharp container



### Step 2:

Continue to fill container until sharps reach the container fill level



### Step 3:

When contents reach the fill line, close and lock for permanent disposal



## Empowering Clinicians

Secure a Sharp products are designed to not only meet the regulatory needs of healthcare, but also complement the patient-care pressures of the clinical environment. Easily relocatable and designed with an inbuilt secure lid closure, **Secure a Sharp containers allow for the safe disposal of sharps in any healthcare setting.**