

How To Use Your Secure a Sharp



SAFELY LOCATE YOUR SECURE A SHARP CONTAINER:



CORRECT DISPOSAL

To ensure your safety in the disposal of medical sharps:

- Do not remove or recap needles
- Do not force sharps into, or overfill the container
- Do not put fingers inside the container

SASH231004

PLEASE NOTE

Please be mindful of your environment and ensure that sharps containers are positioned in a location that minimizes risk for access and disposal.

866 381 8889

WWW.SECUREASHARP.COM

